



HOUND LABS®

HOUND® CANNABIS BREATHALYZER – Collect + Send

PACKAGING COLLECTED SAMPLES

GUIDELINES FOR COLLECTED SAMPLES

This document outlines guidelines for shipping collected HOUND® CANNABIS BREATHALYZER breath samples to the laboratory. For the purposes of these guidelines, HOUND® CANNABIS BREATHALYZER breath samples are defined as non-infectious materials including, but not limited to, methanol. In addition to these guidelines, all shipments must comply with all applicable local, state, and federal laws governing packing, marking, and labeling.

To access resources and training materials, visit houndlabs.com/collector-portal.

GENERAL PACKAGING GUIDELINES

1. Primary Watertight Receptacle

The HOUND® CANNABIS BREATHALYZER cartridge securely houses the methanol-containing inner receptacle in a strong plastic outer casing. Post-collection, the inner receptacle is sealed and held firmly in place via a locking assembly.

2. Laboratory Sample Bag

Following a collection, the collector inserts the compressed post-collection HOUND® CANNABIS BREATHALYZER donor and air cartridges into the supplied laboratory sample bag (along with appropriate documentation) and seals the bag with the integrated adhesive strip.

NOTE: Integrated into the bag are pads intended to absorb any methanol that may leak from either cartridge during shipment.

3. Sturdy Outer Packaging

After sealing, the collector inserts the sealed laboratory sample bag into standard FedEx Express packaging. As possible, cartridges should be cushioned in the shipment by filling empty space with paper, bubble wrap, or other protective materials.

NOTE: Multiple sealed laboratory sample bags may be placed in the same package.

MARKING AND LABELING

Label the package with the supplied, pre-printed shipping label. Once labeled, release to the appropriate delivery service for transport to the laboratory.

NOTICE: Transportation providers may refuse to accept packages that do not meet FedEx, government, or IATA and ICAO requirements. This document is in no way intended to replace requirements mandated by 49CFR and IATA. This guidance is for informational purposes only.

The HOUND® CANNABIS BREATHALYZER is intended to detect recent cannabis use. It does not measure whether, or how much, a person is impaired. It is intended solely for use in employment, law enforcement, and insurance settings.