




HOUND LABS RETRIEVER™ RELEASE NOTES

Release Date	Tuesday, March 26, 2024
Release Version	1.6.310.prod.127
System Downtime	Up to 1 hour beginning at: 7 pm ET; 6 pm CT; 4 pm PT
Customer Impact Risk	Minimal
Communicated Externally?	Proactively communicate and post on collector portal
Downtime Impact	During the outage, users CANNOT: <ul style="list-style-type: none">• Upload collection records to Retriever™ results portal• Receive onboarding information such as user IDs and passwords• Undergo provisioning or interface testing in manufacturing During the outage, users CAN: <ul style="list-style-type: none">• Perform collections
Release Summary	<ul style="list-style-type: none">• Enhancements to user-interface, notifications, and user log-in.

CHANGE: ENHANCED PASSWORD COMPLEXITY

- **Impact:** Upon updating an existing password or creating a new one, users will be notified that passwords must now contain 10 characters with one upper case, one lower case, one number, and one special character. Allowable special characters: !@#\$%^&*
- **Purpose:** This increases application and data security.


New user welcome email with 'RESET PASSWORD' link:





HOUND LABS®

Existing user email with 'RESET PASSWORD' link:




CHANGE: DENY USE OF LAST FIVE (5) PASSWORDS

- **Impact:** A user will not be able to reuse their last five (5) passwords.
- **Purpose:** This increases application and data security.

CHANGE: 'LATEST NEWS' HEADER

- **Impact:** Users now see a 'latest news' header on the Retriever™ landing page.
- **Purpose:** Provides a communication channel for Hound Labs to share news with users about upcoming releases, outages, holidays, etc.





CHANGE: ORG ADMIN + ADMIN MANAGER ROLES CAN RESET USER PASSWORDS

- **Impact:** Users configured with Org Admin or Admin Manager roles can issue reset passwords for other users within their organizations. See 'reset password' email above.
- **Purpose:** Increases application security and allows users to self-serve if a password is compromised or lost.

Organizations > Hound Labs Example

Hound Labs Example

Import/Export : [+ User](#)

Dashboard Organization Hierarchy Locations Results Work Orders Base Stations People

Active ▾ All Loc... ▾ Users ▾ Users for Hound Labs Example

EMAIL ▾	FIRST NAME	LAST NAME	LOCATION	ROLE	BASE STATION ROLE	CREATED	Password Reset
active.user@houndlabs.com	Active	User		Org Admin		03/20/24 3:11pm	reset password

Page 1 < >